




Makes therapy and home practice more productive by providing real-time audio and visual feedback.

Biofeedback is a powerful tool for improving performance during the training or rehabilitation of a variety of motor skills

Real-time Display & Playback are available throughout the learning process to enhance motivation and refine production.

Simple Screens make anyone feel like a technology pro so the focus of practice is *practice* and not “how does this program work?”



Virtual Voice Trainer Practice Time	
Assignment	Minutes
Warm-up and Stretch - Vowel	0.0
Warm-up and Stretch - Nasal	0.0
Warm-up and Stretch - Lip Trill	0.0
Pressed vs. Relaxed	0.0
Starting the Tone - Vowel Onset - Glottal	0.0
Starting the Tone - Vowel Onset - Aspirated	0.0
Starting the Tone - Vowel Onset - Smooth	0.0
Projecting the Voice - Brightness One	0.0
Projecting the Voice - Brightness Two	0.0
Projecting the Voice - Loudness	0.0
Nasality of the Voice	0.0
Assignment One	0.0
Assignment Two	0.0
Open Practice	0.0
Total Practice	0.0

On-Target & Off-Target pre-recorded voices provide concrete models for goal setting.

User-Friendly Interface a hallmark of all Vocal Innovations products

Time Tracker measures client compliance

System Requirements: Windows® 7 through 10

Visit www.EstillVoice.com/Products to learn more.
To order, contact: [info @EstillVoice.com](mailto:info@EstillVoice.com)



55 Standish Boulevard
Pittsburgh, PA 15228 USA
412-341-0316 x701

Voice Training with Vision!